

Presenter PC:

- 1.1
1. Extend view (it's usually duplicated): Windows+P → Erweitern
2. Test keyboard inside the room
3. Sound system:


- a. Switch on speakers and attenuator
- b. Speakers' wheel on "LINE"
- c. PC volume: 40% & Lautsprecher
- d. Speakers volume: white stripes

4. Folder with all scripts:

\\smb.uni-oldenburg.de\psychologie\$\Neuro\data\projects\alo_alpha_neurofeedback\code\experiment

- 1.2
- MATLAB scripts: open in 2 different MATLAB instances (no tabs in one instance) (as in the picture below)

- Online Analysis (ring_buffer.m)
- Experiment (ptb_demo.m)



- 1.3
- Open Purr Data patch: red_noise.pd
- There will be 2 windows, minimize the bigger one and leave the smaller (picture →)

To activate the audio there are 2 options


1. DSP green toggle
2. → Medien → An

If you want to double-check if the audio works:

→ Medien → Teste Audio → OUTPUT MONITOR options

If the audio sounds awful, a possible solution is:

→ Einstellugen → Ausgabe-Geräte → Lautsprecher (Sound Blaster X)





Start

- 2.1 Prepare participant & EEG cap
- Open mbt, connect amplifier and stream signal at 250 Hz
- 2 options to run simulations without participant
1. Open mbt and run the dummy signal
 2. Stream an existing dataset:
 - Open a third matlab instance with new_streamer.m (same folder as ring_buffer.m and ptb_demo.m)
 - Choose dataset to stream (line 7) and press RUN
 - Wait, you'll see a plot (showing O2 channel of participant's EEG) and the Command Window will let you know when it's running.

- 2.2 Run Online Analysis (ring_buffer.m)
- Command Window will tell you when it's ready to start the experiment control script

- 2.3 Run Experiment Control Script (ptb_demo.m)
- A purple screen with a green frog should appear on participant's screen, followed up by a Welcome message asking participants to wait.
 - Start the LabRecorder App once the MATLAB command window tells you.
 - Select all streams: EEG, AAA, EM, NFB
 - Inform participants that they can start by pressing the [spacebar]
 - During the experiment you don't need to do anything: the command window will keep you updated of participant's progress through the instructions and trials.



- 2.4 Purr Data patch troubleshooting:
- If **“error: audio I/O stuck... closing audio”**
- close and open purr-data
 - this won't affect the recording, you don't need to restart.

Click out of the window before starting so when the participant press the [spacebar] they won't toggle the DSP accidentally.

```
DSP ☒
a1.
pdlua: compiled for pd-0.54 on Jun  5 2025 23:38:51
Using lua version 5.4
pdlua: using JavaScript interface (nw.js)
libdir_loader: added 'mrpeach' to the canvas-local objectclass path
rojo~ x[n] = scale * (alpha * x[n-1] + z[n])
scale 0.100000
alpha 0.700000
rojo~ x[n] = scale * (alpha * x[n-1] + z[n])
scale 0.100000
alpha 0.700000
rojo~ x[n] = scale * (alpha * x[n-1] + z[n])
scale 0.100000
alpha 0.700000
rojo~ x[n] = scale * (alpha * x[n-1] + z[n])
scale 0.100000
alpha 0.700000
rojo~ x[n] = scale * (alpha * x[n-1] + z[n])
scale 0.100000
alpha 0.700000
error: audio I/O stuck... closing audio
```


End

3.1 * Important: the following order is relevant to ensure a proper offline reconstruction

1. Stop streaming signals in mbt (smarting streamer app)

2. The online analysis script (ring_buffer.m) should end **automatically**

- The command window will tell you where the chunksize got saved (line 13)
- The plot below should pop-up
- If not, you need to run the last section manually (line 280)



3. Stop LabRecorder

3.2 The experiment control script (ptb_demo.m) Command Window will show where the variables got saved

- Variables: score.mat, trial.mat, task.mat, thresh.mat and storage.mat
- If it didn't run till the end, you need to run the last section to save those variables (line 2350).
- Move out these variables fast, they will get rewritten if you do another recording.

You need to manually close the experiment control script (ptb_demo.m).

- Close the last screen with [Esc].